

Accessories for RaySafe Xi and RaySafe Solo



CASES



1922014/1922016*
RAYSAFE Xi STORM CASE

Heavy duty, waterproof case with customized interior to fit all currently available RaySafe Xi detectors and the RaySafe flexi stand.



1922035/1922037*
RAYSAFE Xi STORM CASE WITH PC INLAY

Storm case with customized interior with space for the RaySafe Xi system and a dedicated space of 275 x 200 x 50 mm, (10.8 x 7.8 x 2.0 inches) for a small laptop or netbook.



1922033/1922036*
RAYSAFE Xi STORM CASE TOSHIBA

Storm case with customized interior for a complete RaySafe Xi system, the RaySafe FPD holder Toshiba and space for a Toshiba specific phantom.



1922013/1922028**
RAYSAFE Xi STANDARD CASE

The standard RaySafe Xi aluminium case with space for a complete RaySafe Xi system and the RaySafe Xi/Solo flexi stand.



1922027
RAYSAFE SOLO STANDARD CASE

The standard RaySafe Solo case has space for the different RaySafe Solo systems and fits the RaySafe OPG holder. The RaySafe flexi stand cannot fit in this case.


SYSTEM OPTIONS:

- * If replacing the standard aluminium case of the RaySafe Xi or RaySafe Solo to a Storm case at time of purchase.
- ** If replacing the standard RaySafe Solo aluminium case to a RaySafe Xi aluminium case at time of purchase


FILTERS



1901017
HVL FILTER SET MAM/RAD
90x90 mm 99.5 % Al filter set
(1 pc 2mm, 2 pcs 1mm, 2 pcs 0.5 mm,
5 pcs 0.1 mm).



1921038
HVL FILTER SET WITH HIGH PURITY
100x100 mm 99.999% Al filter set (6 pcs
0.1 mm). Individually measured thickness
is noted on each filter with a tolerance
of 5%.



1901018
HVL FILTER SET RAD/DENT
90x90 mm 99.5 % Al filter set
(1 pc 2 mm, 2 pcs 1 mm, 2 pcs 0.5 mm)


HOLDERS FOR RAYSAFE Xi/SOLO DETECTORS



1902014
RAYSAFE Xi/SOLO DETECTOR HOLDER
Detector holder for easy placement of
a detector in a bucky or on an image
intensifier.



1902025
RAYSAFE Xi/SOLO OPG HOLDER
Magnetic holder for the detector incl.
positioning device with fluoro screen,
gafchromic film and metal plate.



5600136
GAFCHROMIC FILM XR-QA2 -10PCS
Extra sheets of gafchromic film to be
used together with RaySafe Xi/Solo OPG
holder (1902025).



1902030
RAYSAFE Xi/SOLO FLEXI STAND
For positioning the RaySafe Xi/Solo R/F,
MAM and CT detectors for various
applications. It comes with top, bottom,
two rods, a detector holder, a flexi
coupling, two rubber squares and a 1/4"
camera adapter.



1902037
RAYSAFE Xi/SOLO CASSETTE HOLDER
24x30 cm cassette holder to position the
detector in a bucky.



1902040
RAYSAFE Xi SECTRA DETECTOR HOLDER
Detector holder for easy placement
when measuring on Sectra
mammography systems. To be used
with RaySafe Xi scanning detector.

HOLDERS FOR RAYSAFE Xi/SOLO DETECTORS




1902042
RAYSAFE Xi GE MAMMO DETECTOR HOLDER
 Detector holder which adapts to the holder provided by GE on mammography machines. Measurement position will be identical to ion-chamber measurements.



1902043
RAYSAFE Xi GE R/F DETECTOR HOLDER
 Detector holder which adapts to the holder provided by GE on radiography machines.



1902060
RAYSAFE VACUUM HOLDER
 Detector holder with suction cup for various placements. Attached with a flexible arm and compatible with parts from 1902030 and 1902014.



1902054
RAYSAFE VASCULAR HOLDER GE
 Detector holder for RaySafe Xi/Solo detector which is adapted to the holders provided by GE on angiography machines. Measurement position will be identical to ion-chamber measurements. Aluminum rods are not included.



1902053
RAYSAFE FPD HOLDER TOSHIBA
 Detector holder which centers the position of the RaySafe R/F detector on Toshiba flat panel detector systems; 8x8, 12x12, 16x12 and Kalare. RaySafe Xi Storm case - Toshiba (1922033) has an adapted inlay to fit this detector holder.

PHANTOM



1921046
HEAD AND BODY CT PMMA KIT IN STORM CASE
 Head and body phantom in frosted PMMA for CTDI measurements. The body phantom has a diameter of 32 cm and 4 holes. The head phantom has a diameter of 16 cm and 5 holes. Delivered in a wheel carried, heavy duty Storm case with interior to also fit a RaySafe Xi R/F and CT detector. PMMA material with density 1.19g/cm³.

COMMUNICATION/CONNECTION



1901046
USB SERIAL ADAPTER
 It converts the communication port from USB to RS232 serial interface.



1902028
RAYSAFE Xi mAs CABLE, 10 M
 Long RaySafe Xi mAs cable



1922010
RAYSAFE Xi/SOLO BLUETOOTH SERIAL ADAPTER
 Adapter that fits to the serial port of the RaySafe Xi/Solo base unit for Wireless communication to PC via Bluetooth.



1922012
RAYSAFE mAs RECTIFIER
 For mA/mAs measurements on equipment with an AC mA waveform.



1922019
CNET USB BLUETOOTH ADAPTER CLASS 1
 CNet CBD-120 (100m). Bluetooth adapter for PC with long range communication capability.



1922039
RAYSAFE Xi POWER SUPPLY
 Power supply to be used with the RaySafe Xi/Solo base units.



1922021
RAYSAFE Xi LIGHT GX ADAPTER FOR PHILIPS I.I.
 A Gx adapter for Philips image intensifiers. Connects to the RaySafe Xi light detector. The detector needs to be initiated with the Gx menu at the time of purchase of the Gx adapter.