

# Biomedical

# **DPM1B**

# **Pneumatic Transducer Tester**



## **Technical Data**

The DPM1B pneumatic transducer tester is designed to measure the positive and negative pressures of medical devices in either liquid or gaseous form, and to generate pressure within the  $\pm$  300 mmHg range to assist in repair and quality control.

# **Key Features**

- Battery operated
- Generates and measures positive or negative pressures
- Operates with gas and liquid
- Troubleshooting with 1 % accuracy

# **Technical Specifications**

**Pressure Measurement** 

Operating Range: ± 300 mmHg Pressure Generation Range: ± 300 mmHg Accuracy:  $\pm$  1 % of reading or  $\pm$  1 mmHg Units of Measure: mmHg

**Environmental Requirements** 

Operating Temperature: 10 °C to + 40 °C  $(+50 \,^{\circ}\text{F to} + 104 \,^{\circ}\text{F})$ 

**General Information** 

Display/Control: 0.5 in LCD with LO BATT

indication

Power: 9 V alkaline battery Weight: 10 oz (260 g)



About Fluke Biomedical
Fluke Biomedical is the world's leading manufacturer
of quality biomedical test and simulation products. In
addition, Fluke Biomedical provides the latest medical
imaging and oncology quality-assurance solutions for
regulatory compliance.

Today, biomedical personnel must meet the increasing regulatory pressures, higher quality standards, and rapid technological growth, while performing their work faster and more efficiently than ever. Fluke Biomedical provides a diverse range of software and hardware tools to meet today's challenges.

Fluke Biomedical Regulatory Commitment
As a medical device manufacturer, we recognize and
follow certain quality standards and certifications when
developing our products. We are ISO 9001 certified and
our products are:

• FDA Compliant
• CE Certified, where required
• NIST Traceable and Calibrated
• UL, CSA, ETL Certified, where required



# Biomedical

### **Ordering Information**

2249779: DPM 1B Pneumatic Transducer Tester

### **Standard Accessories**

2572314: User/service manual 2242666: 4-way stop cock adapter 614487: 9 V alkaline battery

### **Optional Accessories**

2242666: 4-way stop cock adapter

2212678: Acrodisc filter

### Fluke Biomedical.

Better products. More choices. One company.

### Fluke Biomedical

PO Box 9090, Everett, WA 98206-9090 USA

### Fluke Biomedical Europe AS

Vegamot 8, N-7048 Trondheim, Norway

### For more information, contact us:

In the U.S.A. (800) 648-7952 or Fax (425) 446-5629 In Europe/M-East/Africa +47 73954700 or Fax +47 73954701

From other countries +1 (425) 347-6100 or Fax +1 (425) 446-5629

Email: sales@flukebiomedical.com Web access: http://www.flukebiomedical.com

©2006 Fluke Biomedical. All rights reserved. Specifications subject to change without notice. Printed in U.S.A.

LITO47